

USER GUIDE

Video Upload Utility

TABLE OF CONTENTS

1	OVERVIEW OF VIDEO UPLOAD PORTAL					
2	PURPOSE OF DOCUMENT					
3	Reg	ISTRATION AND LOG IN	2			
	3.1	Registration	2			
	3.2	Login	3			
4	Ном	IE PAGE	3			
	4.1	Home Page Menu Options	3			
	4.2	Send Invitation Link to Parent/Carer	4			
	4.3	Patient Management Dashboard	5			
	4.4	Patient Videos	6			
5	Par	ENT/CARER FUNCTIONALITY	8			
	5.1	Invitation eMail	8			
	5.2	Parent Device Web Interface	8			
6	Арр	ENDICES	. 11			
	6.1	Version Control	11			

world class expertise 🔶 local care

1 OVERVIEW OF VIDEO UPLOAD PORTAL

The Video Upload Utility securely uploads videos of most formats to a secure management portal for access by the Royal Free Child Development department.

2 PURPOSE OF DOCUMENT

This user guide describes the functionality available within management portal.

3 REGISTRATION AND LOG IN

How to register, log in and the main menu

3.1 Registration

This section describes registration process

3.1.1	System manager	Dear Nurse One Epilepsy,		
	registers user details and user receives registration email	You have been registered as a Consu the Video Management Service.	Itant/Nurse on	
3.1.2	User forename and surname.	Your details are:		
		Forename: Nurse One	3.1.2	
3.1.3	User email, this is also	Surname: Epilepsy	212	
	accessing the system	Email: person@nhs.net	3.1.3	
	5 ,	Temporary Password, gg10#e1+	3.1.4	
3.1.4	Temporary user password, the system will ask for this to be changed on first logon	Assigned Role(s): - Consultant/Nurse 3.1.5 Assigned Site: 1122-Barnet Hospital Activation Status: Active		
3.1.5	Assigned role, hospital site and status should be 'Active'	Please login to the service using this p have to change your password when y Please click here to login into the Vide	assword. You will ou first login. o Management	
3.1.6	Link to portal, user should bookmark for future reference	Service. 3.1.6	-	
		Thank you.		
		The IGspectrum Support Team.		

3.2 Login

This section describes login page



4 HOME PAGE

This section describes options available from the home page.

				Video Management Service			Andrew Trew 🗸
E Home							
O Hospital Management 🗸	PATIENT LIST						
O Patient Management 🐱	Hospital PIN:	All		•	Patient name:	1	
Patient List	DOB:	DD/MM/YYYY			NHS number:		
O Cons/Nurse Management ~	Video dato:				Mobile No:		
Generate Link	Theo dute.				WODIC NO.		
	Search Reset						
		_					
	PIN	Name	DOB	NHS Number		Mobile No	Action
	100	Testing	09/04/2004	222222222		+91 7008948176	View Video

4.1 Home Page Menu Options

This section describes main options available from the home page

4.1.1	Home – displays this page	
4.1.2	Generate Link – allows a link email to be sent to a new parent/carer	Home 4.1.1
4.1.3	Patient Management – Displays the patient management dashboard	 Generate Link 4.1.2 Patient Management ~4.1.3

4.2 Send Invitation Link to Parent/Carer

Clicking 'Generate Link displays the following page:

4.2.1	Enter the email of the parent/carer to receive the invitation	SHARE URL					
		Email ID:* 4.2.1					
4.2.2	Select the site/department from the list (4.2.3) which requires to see the video	Hospital Pin: * Select 4.2.2					
4.2.3	List of departments, configured as part of system setup.	Share Link 4.2.4					
4.2.4	Click button to send email to parent/carer	100-Test					
4.2.5	QR code that allows mobile to open video upload page on mobiledevice of parent/carer						

4.3 Patient Management Dashboard

The patient dashboard lists patients for whom one or more videos have been uploaded

PATIENT LIST			Search fields f	or filtering			
Hospital PIN:	All		•	Patient name:			
DOB:	DD/MM/YYY	ſY		NHS number:			
Video date:	DD/MM/YYY	ſŶ		Mobile No:			
Search	Reset						
Dashboard Display							
PIN	Name	DOB	NHS Number		Mobile No	Action	
100	Testing	09/04/2004	22222222222		+91 7008948176	View Video	

Click 'View Video' button to view patient videos



4.4 Patient Videos

Clicking on 'View Video' displays the following page:





4.4.1	This question appears when a patient is over 16 years old and is to confirm that the patient consents to the on-going retention of videos of them. Management of the retention of videos is dependent on retention policies of the department.	The patient has given their own explicit consent for their videos to be held on the Royal Free video portal. 4.4.1
4.4.2	Uploaded video with standard controls	4.4.2 VIDEO NUMBER: 36 BOOKMARKS DETAILS 4.4.6 Time Title Description
4.4.3	Button to add bookmark. Can be pressed while video is playing or paused. Bookmark entry box will be displayed.	C006 / 030
4.4.4	Title of bookmark	Bookmark X
4.4.5	Description of bookmark. Click save to add bookmark to list 4.4.6	Please do not mention any patient identification details. Title:* 4.4.4 Lescrotion.* 4.4.5
4.4.6	List of bookmarks. Elapsed playing time of video displayed in the 'Time' column to allow navigation back to the relevant point in the video.	Rend

5 PARENT/CARER FUNCTIONALITY

5.1 Invitation eMail

The Parent will receive the following eMail when the 'Generate Link' function is used:



The QR code can be scanned by a mobile device to open the upload web page on the device. Web Pages displayed on device:

5.2 Parent Device Web Interface

Consent page displayed when QR code is scanned

world class expertise 🔶 local care





6 **APPENDICES**

6.1 Version Control

VER	DATE	Αстіо	COMMENTS	STATUS
1.0	05/01/2021	Initial version	First version of user guide	Released

